Immunitics® C6 Lyme ELISA™ Kit
Cat. No.: DK-E601-096
96 Tests
For *In Vitro* Diagnostic Use

INTENDED USE:

The Immunitics® C6 Lyme ELISA™ Kit is an *in vitro* diagnostic test, intended for use as a diagnostic aid in the detection of IgG and IgM antibodies to *Borrelia Burgdorferi* (*B. burgdorferi*) in human serum. The diagnosis of Lyme disease must be made based on history, signs (such as erythema migrans), symptoms, and other laboratory data, in addition to the presence of antibodies to *B. burgdorferi*. Negative results should not be used to exclude Lyme disease.

SUMMARY:

Lyme disease is a multi-system illness caused by infection with any of the three genospecies of the spirochete *B. burgdorferi sensu lato*^{1,2}, including *B. burgdorferi*, *B. garinii*, and *B. afzelii*. Transmission of the disease occurs through the bite of any of several species of *Ixodes* ticks, which are found in the United States³⁻⁵, Europe⁶, Russia, and other countries in Asia.

As it progresses through various stages, Lyme disease produces an array of clinical symptoms^{2,7}. The first sign of infection (primary stage) is the development of a circular skin rash, termed 'erythema migrans' (EM) at the bite site. This occurs in 60-80 % of patients within a few days or weeks of the initial infection; some patients do not develop or may overlook it. General flu-like symptoms (headache, abdominal pain, and fatigue) frequently accompany or follow this rash. Weeks to months later, the disease evolves into its secondary stage. This stage is usually characterized by musculoskeletal pain or arthritis, neurological abnormalities and/or cardiac complications. Brief arthritic attacks affecting large joints become more persistent with time and may resolve into a chronic condition in late infection (tertiary stage). Other manifestations of late stage infection include Acrodermatitis Chronicum Atrophicans, a skin lesion, and various neurological disorders associated with neuroborreliosis.

Serological testing has been demonstrated to be useful in detecting the antibody response to *B. burgdorferi* as an adjunct to diagnosis. In most Lyme patients, an antibody response to spirochetal antigens can be detected within weeks after infection. The sensitivity and specificity of the response as well as the time course and antibody class vary depending on the antigen.

The antigen used in the Immunitics® C6 Lyme ELISA™ kit is a synthetic peptide (C6 peptide*) derived from the VlsE protein, which has been shown to be both specific and highly immunogenic¹³⁻²⁰. The peptide sequence is conserved and equally antigenic in humans infected with *Borrelia burgdorferi sensu stricto* or with European genospecies including *Borrelia afzelii* and *Borrelia garinii*¹⁶. As the antigen represents a defined sequence within the protein, potential cross-reactivity with unrelated and partially related antigens found in other organisms is greatly reduced.

*Protected by patents US 6,475,492, EP1171605, AU767955, CA2370493.

PRINCIPLE:

The Immunitics® C6 Lyme ELISA™ kit is based on a synthetic peptide antigen (C6 peptide) in microwell ELISA format. In the assay procedure, diluted serum samples are added to and incubated in wells of an antigen-coated microwell plate. Antibodies specific to the C6 peptide in the serum sample are bound by the immobilized antigen and unbound antibodies are removed by wash steps. The bound antibodies are detected by addition of a horseradish peroxidase-conjugated (HRP) goat anti-human IgG/IgM conjugate. After removal of excess conjugate by further wash steps, a chromogenic peroxidase substrate containing tetramethylbenzidine (TMB) is added. A blue-green product is produced in wells where antibodies have been bound to the antigen. The color development reaction is quenched by addition of dilute sulfuric acid, producing a yellow color change, after which optical absorbance at 450 nm, corrected by a 590 – 650 nm background subtraction, is measured in each well using an ELISA microplate reader.

REAGENTS SUPPLIED:

1. Microwell Plate (Part # CB-P005-096). 96 wells, provided in twelve strips each containing 8 wells and stored in a resealable foil pouch with desiccant.
2. Positive Control (0.300 mL) (Part # CB-L018-300). Human source material with Gentamicin and ProClin added as preservatives.
3. Negative Control (0.300 mL) (Part # CB-N023-300). Human source material with Gentamicin and ProClin added as preservatives.
4. Cutoff Calibrator (0.300 mL) (Part # CB-N030-300). Human source material with Gentamicin and ProClin added as preservatives. Used to calibrate the cutoff value of the assay.
5. Conjugate (15 mL) (Part # CB-A028-015). Ready to use goat anti-human IgG/IgM horseradish peroxidase conjugate.

6. Sample Diluent (30 mL) (Part # CC-S004-030). Ready to use dilution buffer containing Gentamicin and ProClin as preservatives.
7. TMB ELISA Substrate (15 mL) (Part # CC-S003-015). Ready to use solution containing hydrogen peroxide and TMB.
8. 10X Wash Buffer Concentrate (60 mL) (Part # CC-B001-060). Contains phosphate-buffered saline and Tween-20 detergent, with ProClin added as a preservative.
9. Stop Solution (15 mL) (Part # CC-S005-015). Ready to use acidic solution.

PRECAUTIONS:

1. The human source materials used in this kit have tested negative by FDA approved methods for antibodies to HIV-1 and HIV-2, Hepatitis C and Hepatitis B surface antigen. However, because no test methods can offer complete assurance that those infectious agents are absent, all controls and test specimens should be handled as if capable of transmitting infectious agents. They should be considered as potentially infectious materials and handled at the Biosafety Level 2 as recommended in the CDC/National Institutes of Health manual "*Biosafety in Microbiological and Biomedical Laboratories*, 5th Edition", 2009.
2. Wear proper personal protection while running the assay. Do not allow reagents to come into contact with skin, eyes or mouth, as irritation may result. Wash thoroughly with water in the event of any contact. Do not mouth pipette reagents.
3. Cap controls and calibrator tightly after use to minimize evaporation.
4. Stop Solution contains sulfuric acid. Waste containing Stop Solution should be brought to neutral pH before disposal. Caution: The addition of sodium hypochlorite (bleach) to solutions containing sulfuric acid will produce toxic chlorine gas.
5. Use fresh pipet tips for pipetting each sample and reagent. Re-use of pipet tips may introduce contamination.
6. Dispose of used assay reagents and samples by proper biohazard procedures.
7. Components of kit lots should be used only in the combination supplied, with the exception of 10X Wash Buffer which may be used interchangeably between lots. Do not use kit components beyond the expiration date on the label.
8. All kit components should be brought to room temperature (19 – 26 °C) prior to starting the assay.
9. TMB ELISA Substrate is sensitive to light and metallic ions. Avoid exposure to light and metallic ions.
10. Use only distilled or deionized water for preparation of buffers in the assay.
11. Wear gloves when handling microwell plate strips.
12. For optimal, reproducible and accurate performance of the test, follow instructions in package insert. Deviations from instructions in this insert may lead to false results.
13. Do not pre-wash wells with Wash Buffer or Sample Diluent.
14. Do not mix bottle caps of TMB ELISA Substrate and Conjugate.

MATERIAL REQUIRED BUT NOT SUPPLIED:

1. Graduated cylinders to measure liquid
2. 1 liter flask for Wash Buffer
3. Pipettors to dispense 10 µL to 200 µL (multichannel pipettor recommended or 10 µL, 100 µL and 200 µL pipettes)
4. Uncoated microwell plate or test tubes for sample dilution
5. Distilled or deionized water
6. Timer (0 – 60 minutes)
7. ELISA reader with 450 nm filter and a 590 – 650 nm reference filter
8. Vacuum aspirator or disposal means for assay solutions (multichannel aspirator recommended)
9. Absorbent towels
10. Centrifuge that is capable of 10,000 rpm

OPTIONAL EQUIPMENT:

1. Automated plate washer

STORAGE AND SHELF LIFE OF REAGENTS:

1. Store unopened kit between 2 – 8 °C. Expiration date is shown on the kit label.
2. Store opened kit components between 2 – 8 °C. Components must be used within 60 days of opening.
3. Unused microplate strips should be stored between 2 – 8 °C sealed within the original foil pouch with desiccant.
4. Kit components including Controls, Calibrator, HRP Conjugate, Substrate, Stop Solution, Sample Diluent, and 10X Wash Buffer Concentrate should be stored between 2 – 8 °C.
5. 1X (working) Wash Buffer solution can be stored at room temperature (19 – 26 °C) for up to 60 days.
6. TMB ELISA Substrate should be stored in light protective bottles that do not contain metallic ions. The shelf life of the substrate is based on storage in the plastic amber bottle that is provided with the kit.
7. Do not use kit components beyond the expiration date.

SPECIMEN COLLECTION:

1. Fresh human serum specimens should be collected and either stored at 2 – 8 °C if testing will take place within 10 days, or frozen.
2. Hemolyzed or lipemic sera, and sera subject to multiple freeze-thaw cycles may yield anomalous results.
3. Centrifuge sera for 1 minute at 10,000 rpm in a microcentrifuge to remove any visible precipitate prior to testing.
4. Do not test samples showing evidence of microbial contamination.

5. The assay is intended for use with human serum. The use of plasma with this assay has not been established.

REAGENT, CONTROL, CALIBRATOR, AND SAMPLE PREPARATION:

1. 1X Wash Buffer: To prepare once, add 60 mL of 10X Wash Buffer Concentrate to 540 mL deionized/distilled water in a 1-liter flask and stir thoroughly. **Note:** On removal of the 10X Wash Buffer Concentrate from refrigeration, undissolved salts may be present. Allow the reagent to reach room temperature and shake the bottle to dissolve the salts.
2. Controls and Calibrator: Centrifuge Controls and Calibrator at $10,000 \times g$ for 10 seconds before using to deposit all liquid at the bottom of the tube. (Prior to dispensing, if Controls are noticed to be excessively cloudy, have evidence of clotting or are difficult to pipette, centrifuge a second time for 1 minute at 10,000 rpm.) Dispense 200 μ L of Sample Diluent into 3 clean test tubes or uncoated microplate wells. Add 10 μ L of Negative Control to one tube or well; add 10 μ L of Positive Control to another tube or well; add 10 μ L of Calibrator to a third tube or well and mix thoroughly.
3. Patient samples: Dispense 200 μ L of Sample Diluent into each of a series of tubes or uncoated microplate wells, sufficient for the number of samples to be tested. Add 10 μ L of each sample to the corresponding tube or well and mix thoroughly.

ASSAY PROCEDURE:

Bring all assay reagents to room temperature (19 – 26°C) before beginning the assay. All steps are performed at room temperature (19 – 26°C).

1. Record the sample identity for each well on the provided record sheet to determine the number of strips necessary to perform the assay. Three wells will be needed for controls and calibrator. One well will be needed for each sample.
2. Remove the microplate frame containing the microplate strips from the foil pouch. Remove unneeded microplate strips from the frame and reseal unused strips in the foil pouch with desiccant. The microplate frame should be retained at the end of the assay to be used with the remaining microplate strips.
3. Add 100 μ L diluted Positive Control to one microwell, 100 μ L of diluted Negative Control to another microwell and 100 μ L of diluted Calibrator to a third microwell.
4. Add 100 μ L of each diluted patient sample to microwells.
5. Incubate for 30 minutes.
6. Aspirate wells. If using manual or semi-automated washing manifold, wash three times as follows. Dispense approximately 300 μ L of 1X Wash Buffer into each well, then aspirate. Refill the wells with approximately 300 μ L of 1X Wash Buffer and aspirate a second time. Refill the wells with 300 μ L of 1X Wash Buffer and aspirate a third time. Make sure that all wells have been aspirated after the third (final) wash step. (If an automated plate washer is used, wash four times with each wash consisting of 300-350 μ L of 1X Wash Buffer.) After the final wash for both manual and automated washing, tap the plate on absorbent towels to remove all residual liquid.
7. Dispense 100 μ L Conjugate into each well.
8. Incubate for 20 minutes.
9. Aspirate wells. Wash the wells with 1X Wash Buffer as in Step 6 above.
10. After the final wash tap the plate on absorbent towels to remove all residual liquid.
11. Dispense 100 μ L Substrate into each well.
12. Incubate for 4 minutes. Please note: Optimal assay performance requires precise timing of the Substrate incubation step; timing should be measured from the addition to the first well.
13. Dispense 100 μ L Stop Solution into each well in the same order as substrate was dispensed in previous step. Tap plate gently to mix contents of wells. Read absorbance values within 5 minutes.
14. Read Absorbance at 450 nm with a reference filter of 650 nm using an ELISA plate reader. If the reader is not equipped with a 650 nm filter, the use of an alternate filter between 590 – 650 nm will provide equivalent results.

QUALITY CONTROL:

1. Control values must be within the following ranges in order for the assay to be considered valid:
2. Negative Control $A_{450-650 \text{ nm}}$ must be <0.18 .
3. Calibrator $A_{450-650 \text{ nm}}$ must be <0.18 .
4. Positive Control $A_{450-650 \text{ nm}}$ must be >1.2 .
5. If any control $A_{450-650 \text{ nm}}$ value is not within the above ranges, the assay should be repeated.

CALCULATIONS:

1. Calculate the assay cutoff value by adding 0.3 to the Calibrator absorbance value.
2. Calculate the Lyme Index value (LI) for each patient sample by dividing the $A_{450-650 \text{ nm}}$ of the sample by the cutoff value.

INTERPRETATION OF RESULTS:

Lyme Index

Interpretation

≤ 0.90

Negative result. No antibody to C6 peptide detected in the present assay. This result does not exclude the possibility of *B. burgdorferi s.l.* infection, and where early Lyme disease is suspected, a second sample should be drawn 2 – 4 weeks later and re-tested.

- 0.91 – 1.09 Equivocal result. The imprecision inherent in any method implies a lower degree of confidence in the interpretation of samples with $A_{450-650\text{ nm}}$ values very close to the calculated cutoff value. For this reason an equivocal category has been designated. Samples yielding equivocal results should be retested. If retesting yields an equivocal result, a new serum specimen should be obtained and tested.
- ≥ 1.10 Positive result. Antibody to C6 peptide detected in the present assay.

The cutoff is determined for each assay run by adding 0.3 to the Calibrator absorbance value. In this way, the cutoff is intended to compensate for run-to-run assay variations, which might otherwise affect sensitivity and specificity. The calibrator has been designed to yield an absorbance value reflecting the background reactivity of normal sera.

LIMITATIONS:

1. A negative result does not exclude the possibility of infection with *B. burgdorferi s.l.*. Patients in early stages of Lyme disease and those who have been treated with antibiotics may not exhibit detectable antibody titers. Patients with clinical history, signs or symptoms suggestive of Lyme disease should be re-tested in 2-4 weeks in the event that the initial test result is negative.
2. A positive result is not definitive evidence of infection with *B. burgdorferi s.l.*. It is possible that other disease conditions may produce artifactual reactivity in the assay.
3. This assay should not be used to screen the general population. The predictive value of the assay is a function of the pre-test probability of Lyme disease in the population tested. Hence, only patients with clinical symptoms of Lyme disease or suspected exposure to *B. burgdorferi s.l.* should be tested.
4. Hemolyzed, lipemic, bilirubinemic or turbid samples may produce artifactual assay results. A fresh sample should be collected for re-testing.
5. Although special training is not needed, optimal performance requires strict adherence to the assay procedure described in the insert. Deviations from the procedure may lead to aberrant results.
6. False positive results may be obtained with sera from patients with diseases other than Lyme disease including syphilis, periodontal disease, rheumatoid arthritis, systemic lupus erythematosus, and other autoimmune or infectious diseases.

PERFORMANCE CHARACTERISTICS:

The performance of the C6 Lyme ELISA™ kit was compared with that of the Two-Tier protocol recommended by the CDC¹² (samples found positive or equivocal on ELISA are retested by IgG and IgM Western Blot to determine the final result) with respect to specificity and sensitivity of detection.

SPECIFICITY:

Serum specimens were obtained from 1,842 normal blood donors, comprising 1,329 sera from individuals residing in regions endemic for Lyme disease (northeastern U.S.) and 513 sera from individuals residing in areas considered non-endemic for Lyme disease (southwestern U.S.). Sera were tested on the C6 Lyme ELISA™ kit and by Two-Tier protocol using a whole cell sonicate ELISA in the first step (Table 1). The estimated clinical specificity of the C6 Lyme ELISA™ kit was calculated as 98.8 %, and is statistically equivalent to that of the Two-Tier protocol for both endemic and non-endemic healthy blood donors ($p > 0.05$).

See Appendix: Table 1. Reactivity in Endemic and Non-Endemic Blood Donors

In addition to specimens from endemic and non-endemic areas in North America, serum specimens were obtained from 107 normal blood donors selected from a European population considered non-endemic for Lyme disease (United Kingdom). Sera were tested solely on the C6 Lyme ELISA™ kit, and the estimated clinical specificity has been calculated as 94.4 %.

SENSITIVITY:

A total of 569 well-characterized serum samples were obtained from North American patients who were diagnosed with Lyme disease based on either (1) documented erythema migrans; (2) culture or PCR positive for *B. burgdorferi*; or (3) Lyme arthritis or other symptoms of disseminated Lyme disease. Positive serology was not a criterion for inclusion in this group.

The C6 Lyme ELISA™ kit was evaluated in comparison with the Two-Tier protocol on the Lyme patient group. Overall, the C6 Lyme ELISA™ kit detected 75.0 % of the Lyme disease patients, compared with 51.5 % found positive by the Two-Tier protocol (Tables 2 and 3). The sensitivity of the C6 Lyme ELISA™ kit vs. Two-Tier protocol for detection of Lyme disease patients with positive culture is shown in Table 3.

See Appendix: Table 2. C6 Lyme ELISA™ kit vs. Two-Tier protocol

See Appendix: Table 3. Sensitivity by Sample Definition Criteria

The sensitivity of the C6 Lyme ELISA™ kit for detection of antibodies to *B. burgdorferi* in Lyme disease patient sera was also evaluated with respect to symptom category. A comparison of sensitivity in the C6 Lyme ELISA™ kit and the Two-Tier protocol is shown in Table 4.

See Appendix: Table 4. Sensitivity by Symptom Category

Lyme disease patient samples stratified by time after onset of illness were tested on the C6 Lyme ELISA™ kit vs. the Two-Tier protocol as shown in Table 5.

See Appendix: Table 5. Sensitivity by Time after Onset

An analysis of reactivity in the subset of IgG and IgM Western Blot-positive Lyme disease patients showed that 99.0 % of patients with positive IgG Western Blot (n= 200) and 95.9 % of patients with positive IgM Western Blot (n= 218) were also positive by C6 Lyme ELISA™ kit (Table 6). Conversely, 39.7% of patients positive on C6 Lyme ELISA™ kit (n= 427) were positive on IgM Western Blot, and 37.7 % (n= 426) were positive on IgG Western Blot.

See Appendix: Table 6. C6 Lyme ELISA™ kit vs. Western Blot Results

PROSPECTIVE STUDY RESULTS:

A prospective study was performed on 1,277 serum samples serially received at two reference laboratories for Lyme disease screening tests. Sera were tested on the C6 Lyme ELISA™ kit and with the Two-Tier protocol (Table 7). Concordant results were found in 98.7 % of the specimens tested.

See Appendix: Table 7. Prospective Samples: C6 Lyme ELISA™ kit vs. Two-Tier Results

CROSS-REACTIVE CONDITIONS:

Sera from 366 individuals with disease conditions other than Lyme disease were tested for cross-reactivity in the C6 Lyme ELISA™ kit. Results for seventeen conditions are provided in Table 8. Two sera were found positive by both C6 Lyme ELISA™ kit and Two-Tier protocol, while one additional serum was positive by Two-Tier protocol only (Table 9). The C6 Lyme ELISA™ kit demonstrated an overall specificity of 99.5 % in all categories.

See Appendix: Table 8. Cross-Reactive Conditions: C6 Lyme ELISA™ kit Results

See Appendix: Table 9. Cross-Reactive Conditions: C6 Lyme ELISA™ kit vs. Two-Tier protocol

REPRODUCIBILITY:

Reproducibility was tested on a panel of 8 specimens, including the kit Positive and Negative Controls and 6 specimens representing 2 negative, 2 weakly reactive and 2 positive sera. Reproducibility was assessed in four analyses: intra-assay, inter-assay, inter-lot, and inter-site. Results from each of these analyses is summarized in the following table.

See Appendix: Table 10. Reproducibility: C6 Lyme ELISA™ kit Results

IMMUNETICS, INC. CONTACT INFORMATION:

Immunitics, Inc.
320 Norwood Park South
Norwood, MA 02062-4659 USA
Tel.: 800-227-4765 or 617-896-9100
Fax: 617-896-9110
Email: info@immunitics.com

Contact Immunitics, Inc. for available translations in other languages.

EUROPEAN AUTHORIZED REPRESENTATIVE:

Obelis s.a
Boulevard, General Wahis 53
1030 Brussels, BELGIUM
Tel: (32) 2.732.59.54
Fax: (32) 2.732.60.03
Email: mail@obelis.net

'C6 Lyme ELISA' is a trademark of Immunitics Inc

'Immunitics' is a registered trademark of Immunitics Inc., a wholly-owned subsidiary of Oxford Immunotec USA Inc., which is a wholly-owned subsidiary of Oxford Immunotec Ltd

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IMMUNETICS, INC CONTACT INFORMATION:

Immunitics, Inc
320 Norwood Park South
Norwood, MA 02062-4659 USA
Тел.: 800-227-4765 или 617-896-9100
Факс: 617-896-9110
Эл. почта: info@immunitics.com

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Immunetics C6 Lyme ELISA™ kit - Appendix

TABLES:

Table 1. Reactivity in Endemic and Non-Endemic Blood Donors

	ENDEMIC	NON-ENDEMIC	TOTAL
Immunetics C6 Lyme ELISA™ kit	98.6 %	99.2 %	98.8 %
Two-tier ELISA/WB kits	99.3 %	99.8 %	99.5 %
n=	1329	513	1842

Table 2. C6 Lyme ELISA™ kit vs. Two-Tier protocol

		TWO-TIER ELISA/WB	
n= 569		POS	NEG
C6 Lyme ELISA™	POS	50.5 % 286	24.8 % 141
	NEG	1.2 % 7	23.7 % 135

Table 3. Sensitivity by Sample Definition Criteria

	ALL LYME DISEASE PATIENTS	CULTURE-POSITIVE PATIENTS
Immunetics C6 Lyme ELISA™ kit	75.0 %	65.4 %
Two-tier ELISA + WB kits	51.5 %	39.2 %
n=	569	237

Table 4. Sensitivity by Symptom Category

	ERYTHEMA MIGRANS	ARTHRITIS	FACIAL PALSY	EARLY NEUROLOGICAL	LATE NEUROLOGICAL
Immunetics C6 Lyme ELISA™ kit	66.5 %	98.2 %	85.7 %	88.6 %	100 %
Two-tier ELISA + WB kits	35.2 %	95.6 %	81.0 %	77.3 %	100 %
n=	403	114	21	44	8

Table 5. Sensitivity by Time after Onset

	1-7 Days	8-14 Days	15-21 Days	22-30 Days	> 30 Days
Immunetics C6 Lyme ELISA™ kit	55.4 %	72.9 %	83.7 %	81.6 %	85.1 %
Two-tier ELISA + WB kits	32.0 %	47.1 %	55.1 %	63.2 %	54.1 %
n=	175	70	49	38	181

Table 6. C6 Lyme ELISA™ kit vs. Western Blot Results

		IgG WESTERN BLOT	
		POS	NEG
C6 Lyme ELISA™	n= 525		
	POS	37.7 % 198	43.4 % 228
	NEG	0.4 % 2	18.5 % 97

		IgM WESTERN BLOT	
		POS	NEG
C6 Lyme ELISA™	n= 526		
	POS	39.7 % 209	41.4 % 218
	NEG	1.7 % 9	17.1 % 90

Table 7. Prospective Samples: C6 Lyme ELISA™ kit vs. Two-Tier Results

		TWO-TIER ELISA/WB	
		POS	NEG
C6 Lyme ELISA™	n= 1277		
	POS	3.0 % 38	1.2 % 15
	NEG	0.1 % 1	95.8 % 1223

Table 8. Cross-Reactive Conditions: C6 Lyme ELISA™ kit Results

DISEASE CONDITION	n	POS	NEG
<i>H. pylori</i>	20	0	20
<i>M. pneumoniae</i>	38	0	38
CMV	20	0	20
EBV	20	0	20
HIV	20	0	20
Hepatitis A	25	0	25
Hepatitis B	23	0	23
Hepatitis C	15	0	15
Flu Vaccinee	25	0	25
ANA	20	0	20
RPR	20	1	19
Lipemic	20	0	20
Icteric	20	0	20
Hemolyzed	20	1	19
SLE	20	0	20
RA	20	0	20
RF	20	0	20
Total	366	2	364

Table 9. Cross-Reactive Conditions: C6 Lyme ELISA™ kit vs. Two-Tier protocol

		TWO-TIER ELISA/WB	
n= 366		POS	NEG
C6 Lyme ELISA™	POS	0.5 %	0.0 %
		2	0
	NEG	0.3 %	99.2 %
		1	363

Total Cross-Reactive Conditions

C6: 99.5% negative
 Two-tier: 99.2% negative

Table 10. Reproducibility: C6 Lyme ELISA™ kit Results

SPECIMEN:	Intra-Assay		Inter-Assay CV%	Inter-Lot CV%	Inter-Site CV%
	AVG LI	CV%			
Positive Control	4.74	8.88	11.57	11.37	21.64
Negative Control	0.16	6.59	6.89	16.92	38.73
1	3.36	10.89	16.85	14.05	21.3
2	1.90	8.57	6.02	19.33	20.6
3	4.60	7.66	7.79	14.97	14.94
4	1.31	16.06	14.69	15.49	24.88
5	0.16	10.3	9.77	28.14	32.25
6	0.03	11.72	18.8	52.58	51.53

REFERENCES:

1. Burgdorfer W, Barbour AG, Hayes SF, et al. 1982. Lyme disease- a tick borne spirochetosis? Science. 216:1317-19.
2. Steere AC. 1989. Lyme Disease. N. Engl. J. Med. 321:586-96.
3. Steere AC, Malawista SE. 1979. Cases of Lyme disease in the United States: locations correlated with distribution of *Ixodes dammini*. Ann. Intern. Med. 91:730-3.
4. Levine JF, Wilson ML, Spielman A. 1985. Mice as reservoirs of the Lyme disease spirochete. Am. J. Trop. Med. Hyg. 34:355-60.
5. Burgdorfer W, Lane RS, Barbour AG, Gresbrink RA, Anderson JR. 1985. The western black-legged tick, *Ixodes pacificus*: a vector of *Borrelia burgdorferi*. Am. J. Trop. Med. Hyg. 34:925-30.
6. Barbour AG, Hayes SF. 1986. Biology of *Borrelia* species. Microbiol. Rev. 50(4):381-400.
7. Dattwyler RJ. 1990. Lyme borreliosis: an overview of clinical manifestations. Lab. Med. 21:290-292.
8. Dattwyler, RJ. 1989. Immunodiagnosis of Lyme borreliosis. Rheum. Dis. Clin. N. Am. 15:727-34.
9. Szczepanski A, Benach J. 1991. Lyme borreliosis: host responses to *Borrelia burgdorferi*. Microbiol. Rev. 55: 21-34.
10. Craft JE, Grodzicki RL, Steere AC. 1984. Antibody response in Lyme disease; evaluation of diagnostic tests. J. Inf. Dis. 149: 789-95.
11. Craft JE, Fischer DK, Shimamoto GT, Steere AC. 1986. Antigens of *Borrelia burgdorferi* recognized during Lyme disease. Appearance of a new immunoglobulin M response and expansion of the immunoglobulin G late in the illness. J. Clin. Invest. 72: 504-515.
12. CDC (1995). Recommendations for test performance and interpretation from the Second National Conference on Serologic Diagnosis of Lyme Disease. Mort. Morbid. Weekly Rep. 44: 590-1.
13. Liang, F.T., Alvarez, A., Gu, Y., Nowling, J.M., Ramamoorthy, R., Philipp, M.T. 1999. An Immunodominant conserved region within the variable domain of VlsE, the variable surface antigen of *Borrelia burgdorferi*. J. Immunol. 163: 5566-5573.
14. Liang, F.T., Steere, A.C., Marques, A.R., Johnson, B.J.B., Miller, J.N., Philipp, M.T. 1999. Sensitive and specific serodiagnosis of Lyme disease by enzyme-linked immunosorbent assay with a peptide based on an immunodominant conserved region of *Borrelia burgdorferi* VlsE. J. Clin. Microbiol. 37(12): 3990-3996.
15. Liang, F.T. and Philipp, M.T. 1999. Analysis of antibody response to invariable regions of VlsE, the variable surface antigen of *Borrelia burgdorferi*. Infection and Immunity 67(12): 6702-6706.
16. Liang, F.T., Aberer, E., Cinco, M., Gern, L., Hu, C.M., Lobet, Y.N., Ruscio, M., Voet, P.E., Jr., Weynants, V.E., Philipp, M.T. 2000. Antigenic conservation of an immunodominant invariable region of the VlsE lipoprotein among European pathogenic genospecies of *Borrelia burgdorferi* SI. J. Infectious Diseases 182.
17. Mogilyansky E, Loa CC, Adelson ME, Mordechai E, Tilton RC. 2004. Comparison of Western immunoblotting and the C6 Lyme antibody test for laboratory detection of Lyme disease. Clin. Diagn. Lab. Immunol. 11(5):924-9.
18. Philipp MT, Wormser GP, Marques AR, Bittker S, Martin DS, Nowakowski J, Dally LG. 2005. Clin. Diagn. Lab. Immunol. 12(9):1069-74.
19. Philipp MT, Marques AR, Fawcett PT, Dally LG, Martin DS. C6 test as an indicator of therapy outcome for patients with localized or disseminated lyme borreliosis. 2003. J Clin Microbiol. 41(11):4955-60.
20. Bacon RM, Biggerstaff BJ, Schriefer ME, Gilmore RD Jr, Philipp MT, Steere AC, Wormser GP, Marques AR, Johnson BJ. 2003. Serodiagnosis of Lyme disease by kinetic enzyme-linked immunosorbent assay using recombinant VlsE1 or peptide antigens of *Borrelia burgdorferi* compared with 2-tiered testing using whole-cell lysates. J Infect Dis. 187(8):1187-99.
21. Wormser GP, Schriefer M, Aguero-Rosenfeld ME, Levin A, Steere AC, Nadelman RB, Nowakowski J, Marques AR, Johnson BJ, Dumler JS. 2013. Single-tier testing with the C6 peptide ELISA kit compared with two-tier testing for Lyme disease. Diagn Microbiol Infect Dis. 75(1):9-15.